

Anti-VCP antibody [Hs-14] (STJ16100810)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/ICC
Host / Source	Mouse
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	7415
Gene Symbol	VCP
UniProt ID	TERA_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---