

Anti-Alpha-Fodrin antibody [AA6] (STJ16100815)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB/ELISA
Host / Source	Mouse
Reactivity	Human/Rat/Chicken/Rabbit/Guinea Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-F WB 1:1000 ELISA
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---