

Anti-Erythroid Alpha-Spectrin antibody [AF10] (STJ16100820)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB/ICC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	6708
Gene Symbol	SPTA1
UniProt ID	SPTA1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------