

Anti-CD9 antibody [MM2/57] (STJ16100822)

GENERAL INFORMATION

| | |
|---------------|---|
| Product Type | Primary antibodies |
| Applications | WB/FC/IHC-F |
| Host / Source | Mouse |
| Reactivity | Human/Llama/Mink/Mustelid/Ferret/Dog/Horse/Bovine/Cat/Rhesus/Pig/Rabbit |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1.0 mg/mL |
| Conjugation | Unconjugated |
| Formulation | PBS with 0.09% sodium azide |
| Isotype | IgG2b |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

TARGET INFORMATION

| | |
|-------------|---------------------------|
| Gene ID | 928 |
| Gene Symbol | CD9 |
| UniProt ID | CD9_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|