

Anti-Monocyte Antigen antibody [IPO-M6] (STJ16100825)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (0.5-1.0 µg/million cells in 0.1 ml). Immunohistology (1-2 µg/ml for 30 min at RT). Immunoprecipitation (1-2 µg/million cells in 1ml anti mouse coated Sepharose-4B suspension).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

