

## Anti-HCAM antibody [DF 1485] (STJ16100826)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/FC/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-F IHC-P 1:5-1:10 FC WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---