

## Anti-Golgi complex antibody [AE-6] (STJ16100835)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/IF/IHC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (1-2 µg/million cells in 0, 1 ml, fix cells in 4% PFA for 10 min, at 4°C, permeabilized with 0, 2% saponin or digitonin for 15 min, at 4°C). Immunoblotting (1-2 µg/ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (1-2 µg/ml f
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1 Kappa
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">10466</a>
Gene Symbol	<a href="#">COG5</a>
UniProt ID	<a href="#">COG5_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

