

Anti-Poly Polymerase antibody [A6.4.12] (STJ16100843)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IP/IHC-F/ELISA/IF
Host / Source	Mouse
Reactivity	Human/D.melanogaster/Hamster/Mouse/Rat/Xenopus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---