

Anti-Actin Alpha-muscle antibody [HHF35] (STJ16100854)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	EM/ELISA/ICC/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Chicken/Monkey/Rabbit/Rat/Swine/Zebrafish

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

