

Anti-Prolactin-epitope C antibody [INN-hPRL-1] (STJ16100865)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-P/IHC-F/ELISA |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1.0 mg/mL |
| Conjugation | Unconjugated |
| Formulation | PBS with 0.09% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|