

Anti-TGF Beta antibody [TB21] (STJ16100866)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA/FC/IHC-F
Host / Source	Mouse
Reactivity	Human/Mink/Sheep/Rabbit/Rat/Mustelid/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Formulation	PBS
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---