

Anti-SKP2 antibody [SKP2-8D9] (STJ16100869)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IF/ELISA/IP
Host / Source	Mouse
Reactivity	Human/Non-Human Primates

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Dilution Range	WB IHC ICC IF 1:250 ELISA IP
Formulation	PBS with 0.1% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	6502
Gene Symbol	SKP2
UniProt ID	SKP2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---