

Anti-Ovarium Carcinoma Antigen antibody [OV632] (STJ16100873)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 µg/mL
Conjugation	Unconjugated
Dilution Range	IHC-F 1:50 FC
Formulation	PBS with 0.01M Sodium Azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

