

Anti-IgG1 Negative Control antibody [CT6] (STJ16100881)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Guinea Pig/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein G purified
Formulation	PBS with 0.1% BSA and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---