

## Anti-MRP1 antibody [MRPm6] (STJ16100887)

### GENERAL INFORMATION

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | FC/WB/ICC/IHC-F/IHC-P |
| Host / Source | Mouse                 |
| Reactivity    | Human                 |

### PRODUCT PROPERTIES

|                     |                                                                          |
|---------------------|--------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                               |
| Concentration       | 250 Mu g/mL                                                              |
| Conjugation         | Unconjugated                                                             |
| Dilution Range      | FC<br>WB<br>ICC<br>IHC-F<br>IHC-P 1:20-1:50                              |
| Formulation         | Serum free tissue culture supernatant with 07% BSA and 0.1% sodium azide |
| Isotype             | IgG1                                                                     |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.                           |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">4437</a>       |
| Gene Symbol | <a href="#">MSH3</a>       |
| UniProt ID  | <a href="#">MSH3_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

