

Anti-MRP1 antibody [MRPr1] (STJ16100888)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/ICC/IHC-F/IHC-P
Host / Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 Mu g/mL
Conjugation	Unconjugated
Dilution Range	FC WB ICC IHC-F IHC-P 1:20-1:50
Formulation	Serum free tissue culture supernatant with 07% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4437
Gene Symbol	MSH3
UniProt ID	MSH3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---