

Anti-MRP1 antibody [MRPm5] (STJ16100889)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/ICC/FC/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Purification	Protein G purified
Dilution Range	IHC-F IHC-P 1:20 ICC FC WB 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4437
Gene Symbol	MSH3
UniProt ID	MSH3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---