

## Anti-MDR antibody [M2I-4] (STJ16100891)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/FC/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 Mu g/mL
Conjugation	Unconjugated
Dilution Range	ICC IHC-F 1:20 FC WB 1:20
Formulation	Serum free tissue culture supernatant with 07% BSA and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------