

Anti-LRP/MVP antibody [LMR5] (STJ16100893)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 Mu g/mL
Conjugation	Unconjugated
Formulation	supernatant with 1% BSA and 0.1% sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	9961
Gene Symbol	MVP
UniProt ID	MVP_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---