

Anti-MRP3 antibody [M3II-21] (STJ16100897)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 µg/mL
Conjugation	Unconjugated
Dilution Range	ICC 1:20-1:50 IHC 1:20 WB 1:20-1:50
Formulation	supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	8714
Gene Symbol	ABCC3
UniProt ID	MRP3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---