

Anti-MRP2 antibody [M2II-12] (STJ16100898)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/FC/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 Mu g/mL
Conjugation	Unconjugated
Dilution Range	ICC 1:20-1:50 IHC-F 1:20 FC WB 1:50
Formulation	supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	1244
Gene Symbol	ABCC2
UniProt ID	MRP2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---