

Anti-Breast Cancer Resistance Protein antibody [BXP-34] (STJ16100903)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-F 1:10-1:20 FC
Formulation	PBS with 1% BSA and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	9429
Gene Symbol	ABCG2
UniProt ID	ABCG2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---