

## Anti-MVP antibody [1011] (STJ16100909)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IF/IHC/IP       |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                                                                            |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                                                                 |
| Concentration       | Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.                                                                                                                                                                                                         |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                               |
| Purification        | Affinity purified from tissue culture.                                                                                                                                                                                                                     |
| Dilution Range      | Immunoblot (1-2 µg/ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (1-2 µg/ml for 30-60 minutes at RT; information on a suitable antigen retrieval method for staining of formalin-fixed tissues is unavailable to date). Immunoprecipitation ( |
| Formulation         | PBS with 0.02% Sodium Azide.                                                                                                                                                                                                                               |
| Isotype             | IgG1k                                                                                                                                                                                                                                                      |
| Storage Instruction | Store for up to 1 year at 2-8°C upon receipt.                                                                                                                                                                                                              |

### TARGET INFORMATION

|             |                           |
|-------------|---------------------------|
| Gene ID     | <a href="#">9961</a>      |
| Gene Symbol | <a href="#">MVP</a>       |
| UniProt ID  | <a href="#">MVP_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

