

Anti-MRP6 antibody [M6II-7] (STJ16100913)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/FC/WB
Host / Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 Mu g/mL
Conjugation	Unconjugated
Dilution Range	ICC IHC-F 1:20-1:50 FC WB 1:20-1:50
Formulation	supernatant with 1% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	368
Gene Symbol	ABCC6
UniProt ID	MRP6_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---