

Anti-Breast Cancer Resistance Protein antibody [BXP-9] (STJ16100916)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB/FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 Mu g/mL
Conjugation	Unconjugated
Dilution Range	ICC 1:20-1:50 IHC-F 1:20 WB 1:50 FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	9429
Gene Symbol	ABCG2
UniProt ID	ABCG2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---