

Anti-Breast Cancer Resistance Protein antibody [BXP-53] (STJ16100917)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB/FC
Host / Source	Rat
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Purification	Protein G purified
Dilution Range	ICC 1:20-1:50 IHC 1:20 WB 1:20-1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	9429
Gene Symbol	ABCG2
UniProt ID	ABCG2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---