

Anti-CD13 antibody [WM15] (STJ16100918)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/IHC-F
Host / Source	Mouse
Reactivity	Human/Non-Human Primates

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	290
Gene Symbol	ANPEP
UniProt ID	AMPN_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---