

Anti-Endothelial Cell Protein Receptor antibody [LMR-42] (STJ16100919)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC
Host / Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-F 1:20 FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---