

Anti-MRP2 antibody [M2III-5] (STJ16100920)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Rat/Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 µg/mL
Conjugation	Unconjugated
Dilution Range	ICC IHC-F IHC-P 1:20 WB
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	1244
Gene Symbol	ABCC2
UniProt ID	MRP2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---