

Anti-MRP4 antibody [M4I-10] (STJ16100921)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB/FC
Host / Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Purification	Protein G purified
Dilution Range	ICC 1:20-1:50 IHC 1:20 WB 1:20-1:50
Formulation	Tissue culture supernatant with 0.1% BSA and 0.09% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	10257
Gene Symbol	ABCC4
UniProt ID	MRP4_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---