

## Anti-TNF-RI antibody [HM104] (STJ16100930)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-F/FC/ELISA/IP  |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

### PRODUCT PROPERTIES

|                     |                                                |
|---------------------|------------------------------------------------|
| Clonality           | Monoclonal                                     |
| Concentration       | 100 Mu g/mL                                    |
| Conjugation         | Unconjugated                                   |
| Formulation         | PBS with 0.1% BSA and 0.02% sodium azide       |
| Isotype             | IgG2a                                          |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">7132</a>        |
| Gene Symbol | <a href="#">TNFRSF1A</a>    |
| UniProt ID  | <a href="#">TNR1A_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|