

Anti-MRP4 antibody [M4I-80] (STJ16100936)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ICC/IHC-F/WB/FC |
| Host / Source | Rat |
| Reactivity | Human/Mouse/Monkey |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 200 Mu g/mL |
| Conjugation | Unconjugated |
| Dilution Range | ICC IHC-F 1:20 WB 1:50 FC |
| Formulation | cell culture supernatant with 0.7% BSA and 0.1% sodium azide |
| Isotype | IgG2a |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 10257 |
| Gene Symbol | ABCC4 |
| UniProt ID | MRP4_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|