

## Anti-Leukocyte Common Antigen antibody [BRA-55] (STJ16100941)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/ELISA/IF/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-F IHC-P 1:10-1:20 ELISA IF WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">5788</a>
Gene Symbol	<a href="#">PTPRC</a>
UniProt ID	<a href="#">PTPRC_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---