

Anti-hRPP20 Immunoglobulin antibody [14B11] (STJ16100943)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB/IF/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg BSA
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---