

Anti-PTPRR protein isoforms antibody [6A6] (STJ16100949)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB/IP
Host / Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg BSA
Isotype	IgG2
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	5801
Gene Symbol	PTPRR
UniProt ID	PTPRR_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---