

## Anti-Citrullinated Fibrinogen Immunoglobulin antibody [3D1] (STJ16100950)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB/ELISA
Host / Source	Mouse
Reactivity	Mouse/Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	ELISA 1:2000-1:4000 IHC 1:100-1:500 WB 1:100-1:500
Formulation	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---