

## Anti-ICAM-1 antibody [MEM-111] (STJ16100956)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human/Cow/Rat      |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Conjugation         | Unconjugated   |
| Formulation         | Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide |
| Isotype             | IgG2a  |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.                                     |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|