

Anti-Mucin-1 antibody [BC-2] (STJ16100957)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-F IHC-P 1:20 WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4582
Gene Symbol	MUC1
UniProt ID	MUC1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---