

Anti-Prostate Specific Antigen antibody [ER-PR8] (STJ16100961)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/IF/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-F IHC-P 1:10 IF WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	354
Gene Symbol	KLK3
UniProt ID	KLK3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---