

Anti-MRP5 antibody [M5I-10] (STJ16100966)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB/FC
Host / Source	Rat
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 µg/mL
Conjugation	Unconjugated
Dilution Range	ICC IHC-F 1:20 WB 1:50 FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	10057
Gene Symbol	ABCC5
UniProt ID	MRP5_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---