

Anti-MRP7 antibody [M7I-3] (STJ16100967)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB/FC
Host / Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	250 µg/mL
Conjugation	Unconjugated
Dilution Range	ICC IHC-F 1:20 WB FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	89845
Gene Symbol	ABCC10
UniProt ID	MRP7_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------