

## Anti-P53 antibody [BP53+12] (STJ16100968)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/IF
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC 1:5-1:10
Formulation	Contains 0.09% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">7157</a>
Gene Symbol	<a href="#">TP53</a>
UniProt ID	<a href="#">P53_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------