

Anti-Lymphocytes antibody [CT-4] (STJ16100971)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC
Host / Source	Mouse
Reactivity	Guinea Pig/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	500 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-F 1:100-1:1000 FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG3
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---