

Anti-MHC class II antibody [CT-10] (STJ16100974)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Guinea Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	700 ul/mL
Conjugation	Unconjugated
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---