

## Anti-MRP6 antibody [M6II-68] (STJ16100976)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	50 Mu g/mL
Conjugation	Unconjugated
Dilution Range	ICC IHC-F 1:20-1:50 WB
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">368</a>
Gene Symbol	<a href="#">ABCC6</a>
UniProt ID	<a href="#">MRP6_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---