

Anti-Fibronectin antibody [568] (STJ16100986)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/IF/IHC
Host / Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (1-2 µg/million cells in 0, 1 ml). Immunoblotting (1 µg/ml for 2h at RT). Immunofluorescence (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30-60 min at RT; staining of formalin-fixed tissues requires boiling tissue sections in 10mM Tris
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2335
Gene Symbol	FN1
UniProt ID	FINC_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

