

## Anti-Human IgA antibody [EBS-huA] (STJ16101010)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/WB/IHC/IP
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT). Flow cytometry (0.5-1.0 µg/million cells in 0.1 ml). Immunoblotting (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; staining of formalin-fixed tissues is served
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">1308</a>
Gene Symbol	<a href="#">COL17A1</a>
UniProt ID	<a href="#">COHA1_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

