

Anti-Kappa Light Chain antibody [EBS-huKappa] (STJ16101014)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: not known; tracer: 0, 001-100 µg/ml for 30 min at RT). Flow Cytometry (0.5-1.0 µg/million cells in 0.1 ml). Immunofluorescence (0.5-1 µg/ml). Immunohistology (formalin-fixed: 1-2 µg/ml for 30 min at RT; requires boiling tissue sec
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene Symbol	IGKC
UniProt ID	IGKC_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

