

Anti-Campylobacter jejuni antibody [EBS-I-002] (STJ16101047)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/IHC
Host / Source	Mouse
Reactivity	C.Jejuni/C.coli/C.hyointestinalis/C.lardis/S.aureus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: not known; tracer: 0.001-100 µg/ml for 30 min at RT). Immunofluorescence (0.5-1.0 µg/ml). Immunohistology (1.0-2.0 µg/ml for 30 min at RT; no suitable antigen retrieval method for formalin fixed material is known to date).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgMk
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---