

Anti-EBV-LMP1 antibody [EBS-I-025] (STJ16101053)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/WB/IP/IHC
Host / Source	Mouse
Reactivity	EBV

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT). Immunofluorescence (0.5-1 µg/ml). Immunoblotting (1.0-2.0 µg/ml). Immunoprecipitation (1-2 µg per 100-500 µg cell lysate protein/1 ml of anti-mouse coated Sepharose-4B su
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

